



Catalogue-2024

Pain Management Aids
Rehabilitation Aids
Immobilizers
Traction Kits



www.tescaglobal.com



Our Family stands behind Tesca Technologies
so you can stand up for yours



Our mission is to make high quality orthopaedic support braces, fracture aids and pain management aids easily accessible because we believe that everyone deserves good health.

We strictly follow international standards and source only the highest quality materials to develop our products.

Our commitment to supporting good health for everyone drives every facet of our business and partnerships. Innovation, dedication and strong relationships with our distribution partners are the core values that define Tesca Technologies Pvt. Ltd.

05

ABDOMEN & SPINE SUPPORT

Abdominal Belt
Tummy Trimmer
Maternity Belt
Ash Brace (Hyper Extension Brace)
Rib Belt
Chest Binder
Lumbo Sacral Support
Contoured L.S. Support
Dorso Lumbar Brace
Hernia Belt
Posture Corrector
Scrotal Support

08

CERVICAL SUPPORT

Cervical Collar Hard
Cervical Collar Soft
Cervical Collar Soft With Support
Cervical Orthosis
(Philadelphia Collar)
Cervical Pillow (Memory Foam)
Cervical Pillow Regular

10

SHOULDER & CLAVICLE SUPPORT

Elastic Shoulder Immobilizer
Universal Shoulder Immobilizer
Clavicle Brace
Clavicle Brace With Buckle

11

ARM & ELBOW SUPPORT

Pouch Arm Sling (Tropical)
Pouch Arm Sling (Baggy)
Elbow Support
Tennis Elbow Support
Arm Immobilizer **NEW/ARRIVAL**
Forearm Splint

14

WRIST & FOREARM SUPPORT

Wrist Brace
Wrist Brace With Thumb Support
Wrist & Forearm Splint (Right/Left)
Elastic Wrist Splint (Right/Left)

15

FINGER SPLINTS

Finger Cot
Finger Extension Splint
Frog Splint
Mallet Finger Splint
Thumb Spica Splint
Garter Splint **NEW/ARRIVAL**

17

TRACTION KITS

Pelvic Binder
Pelvic Traction Kit (with Weight Beg)
Traction Weight Bag
Leg Traction Brace
Cervical Traction Kit + Weight Bag
(Sleeping / Sitting)
Skin Traction Set (PUF Liner)

19

THIGH, KNEE, ANKLE & CALF SUPPORT



- Thigh Support
- Knee Cap
- Knee Cap Open Patella
- Knee Cap Rigid Hinge
- Elastic Knee Support
- Patellar Support
- Knee Immobilizer (14")
- Knee Immobilizer (19")
- Knee Immobilizer (22")
- R.O.M. Knee Brace (22")
- Calf Support
- Foot Drop Splint (Right/Left)
- Ankle Brace
- Ankle Binder
- Anklet

23

COMPRESSION GARMENTS



- Compression Stocking Mid Thigh
- Compression Stocking Below Knee
- Medical Compression Stocking Mid Thigh-Classic
- Medical Compression Stocking Below Knee-Classic
- Prosthetic Compression Stocking Mid Thigh
- Prosthetic Compression Stocking Below Knee

25

FOOT, HEEL & TOES SUPPORT



- Arch Support
- Toe Separator
- Insole
- Heel Cushion
- Metatarsal Pad
- Cast Shoe
- Bunion Splint **NEWARRIVAL**

27

NEOPRENE PRODUCTS



- Shoulder Support (Neoprene)
- Waist Belt (Neoprene)
- Wrist Brace (Neoprene)
- Wrist Brace With Thumb Support (Neoprene)
- Knee Wrap (Neoprene)
- Knee Wrap Hinged (Neoprene)

29

PHYSIO AIDS



- Exercise Band (Actiband)
- Exercise Gel Ball (Soft/ Medium/Hard)
- Weight Cuff

30

WELLNESS PRODUCTS



- Back Rest
- Back Rest - Lumbo **NEWARRIVAL**
- Coccyx Cushion Seat
- Orthopaedic Donut Pillow w **NEWARRIVAL**
- Travel Neck Pillow **NEWARRIVAL**
- Tourniquet
- Orthopaedic Heating Pad **NEWARRIVAL**

32

HOT & COLD THERAPY AIDS



- Hot & Cold Pack

33

OTHER PRODUCTS



- Acticrepe
- (Cotton Crepe Bandage B.P.)
- Disposable 3 Ply Face Mask
- Actipop
- (Plaster of Paris Bandage)



ABDOMEN & SPINE SUPPORT



TD-101


Abdominal Belt

Our abdominal belt is perfect for postoperative conditions of the abdomen and post-delivery umbilical And ventral hernias.

Indications:

- Postoperative conditions of the abdomen.
- Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the navel



TD-107


Tummy Trimmer

Optimal compression of our Tummy Trimmer binds surgicals incisions and helps heal abdomen and chest.

Indications:

- Postural deformity.
- Osteoporosis of the spine.
- Postoperative care of the spine.
- Mild thoraces and lumbar injuries.
- Degenerative Disc Disease
- Herniated Disc or PIVD.

Size	Inches	Cm
S	28-32	70-80
M	32-36	80-90
L	36-40	90-100
XL	40-44	100-110
XXL	44-48	110-120
XXXL	48-52	120-130
XXXXL	52-56	130-140



*Measure circumference around the navel



TD-111


Maternity Belt

This maternity belt provides relief of hip and Pelvic pain through sufficient stretching and gentle compression to reduce bladder pressure and spread the baby's weight.

Indications:

- Postoperative conditions of the abdomen.
- Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
S	32-38	80-95
M	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3" below the navel



TD-108

Ash Brace (Hyper Extension Brace)

Our ash brace is aerodynamically shaped with telescoping, aluminium section, and malleable horizontal main body to form a cross bar.

Indications:

- Osteoporotic Bone Collapse.
- Ankylosing Spondylitis.
- Conservative management of Vertebral Bone Tuberculosis.
- Postoperative, post-traumatic, or chronic spinal pathology.
- Complications of thoracolumbar junction which needs moderate Immobilization.

Size	Feet	Cm
Short	5-6	152-183
Long	Above 6	>183



*Measure height of the patient

ABDOMEN & SPINE SUPPORT



TD-104


Rib Belt

Designed with a splinting pad in comfort, our rib belt is useful in aiding the healing process in rib fracture, relieving chest pain, and intercostal muscle strain.

Indications:

- Postoperative conditions of the abdomen.
- Post-delivery.
- Umbilical and ventral hernia.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the chest below underarms



TD-105


Chest Binder

Designed for comfort with an easy hook & loop closure, moisture-wicking fabric and splinting pad.


Indications:

- Post cardiothoracic surgeries.
- Fracture of the sternum.
- Pain due to chest injury.
- Rib fractures.
- Osteoporosis.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the chest below underarms



TD-102

Lumbo Sacral Support

Our lumbo-sacral support is used to comfortably maintain correct posture and support the lower back while sitting, walking or driving.

Indications:

- Postoperative conditions of the abdomen.
- Post-delivery umbilical and ventral hernia.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the navel



TD-103


Contoured L.S. Support

This Contoured Lumbo-Sacral Support is ideal for postoperative care of lumbar spine and for lower back pain, strain or sprains.

Indications:

- Postoperative care of the lumbar spine.
- Low back pain.
- Low back sprains/strains.
- Lumbar disc injury.
- Degenerative Disc Disease.
- Herniated Disc or PIVD.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the navel

ABDOMEN & SPINE SUPPORT



TD-106


Dorso Lumbar Brace

Features anatomically contoured splints for precise immobilization and clavicle straps with cotton stockinet to correct shoulder posture alignment.

Indications:

- Postural deformity.
- Osteoporosis of the spine.
- Postoperative spine care.
- Mild lumbar and thoraces injury.
- Degenerative Disc Disease.
- Herniated Disc or PIVD.

Size	Inches	Cm
Short	4.6-5.6Ft.	140-171
Long	5.6-6.6Ft.	171-201



*Measure circumference around the navel

Size	Inches	Cm
Universal	28-44	71-112
Special	44-52	112-132



*Measure circumference around the navel



TD-109


Hernia Belt

Designed for a comfortable and supportive fit our hernia belt is lightweight and features adjustable pelvic and leg straps.

Indications:

- Reducible inguinal or scrotal hernia.
- Post hernia repair surgery.
- Post-delivery umbilical and ventral hernia.
- Post abdominal hernia repair surgery.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130



*Measure circumference 3" below navel



TD-110


Posture Corrector

Designed to correct slouching and improve posture.

Indications:

- Back pain.
- Spine misalignment.
- Lumbar pain.
- Poor posture.

Size	Inches	Cm
CH	22-28	56-71
S	28-34	71-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130



*Measure circumference around the chest below underarms



TD-1103


Scrotal Support

Fabricated extra-soft and hypoallergenic materials, our scrotal support is well-ventilated to allow for ease of use and extended wear comfort.

Indications:

- Palliative care.
- Treatment for postoperative care.
- Tenderness, pain or swelling in scrotal area.


Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130



*Measure circumference 3" below navel



CERVICAL SUPPORT



Cervical Collar Hard

Our hard Cervical Collar is made of high-quality polyethylene with padded edges moulded with Ethafoam.

Indications:

- Cervical Spondylitis.
- Cervical Disc Disease.
- Fracture of cervical spine.
- Whiplash injuries.

TD-205




Cervical Collar Soft

Made of soft density foam for gentle and cushioned neck support.

Indications:

- Cervical Trauma.
- Neck Pain and Stiffness.
- Excessive Mobility.
- Rheumatoid Arthritis.
- Strain of Neck Muscles.
- Whiplash Recovery.
- Degenerative Disc Disease.

TD-202




Cervical Collar Soft With Support

Made of firm density foam covered with a soft, smooth cotton stockinet splint inside for better neck support.

Indications:

- Cervical Trauma.
- Neck Pain and stiffness.
- Excessive mobility.
- Rheumatoid Arthritis.
- Strain of neck muscles.
- Whiplash recovery.
- Degenerative disc disease.

TD-201



Cervical Orthosis (Philadelphia Collar)

Made of Ethafoam material to padded support our Philadelphia Collar has a trachea opening in front.

Indications:

- Cervical trauma.
- Cervical sprain.
- Cervical Spondylitis.


TD-204

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
M	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25




*Measure distance between edge of sternum & lower jaw with neck in the extension position

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
M	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50




*Measure distance between edge of sternum & lower jaw with neck in the extension position

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
M	3.5-4.0	8.75-10.00
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50



*Measure distance between edge of sternum & lower jaw with neck in the extension position

Size	Inches	Cm
S	3.0-3.5	7.5-8.75
M	3.5-4.0	8.75-10
L	4.0-4.5	10.00-11.25
XL	4.5-5.0	11.25-12.50



*Measure distance between edge of sternum & lower jaw with neck in the extension position



TD-203

Cervical Pillow (Memory Foam)

Made of memory foam to provide proper support to the neck region while relaxing or sleeping.

Indications:

- Temporomandibular disorders.
- Cervical trauma.
- Morning hand stiffness and swelling.
- Whiplash injuries.
- Arthritis.
- Neck pain.

Size	Inches	Cm	
Universal	3.0-4.5	7.50-11.25	

*Measure distance between edge of sternum & lower jaw with neck in the extension position



TD-206

Cervical Pillow (Regular)

Provides sustained support to cervical vertebrae during sleep.

Indications:

- Temporomandibular disorders.
- Cervical trauma.
- Morning hand stiffness and swelling.
- Whiplash injuries.
- Arthritis.
- Neck pain.

Size	Inches	Cm	
Universal	3.0-4.5	7.50-11.25	

*Measure distance between edge of sternum & lower jaw with neck in the extension position



SHOULDER & CLAVICLE SUPPORT



TD-302


Elastic Shoulder Immobilizer

Our Elastic Shoulder Immobilizer is made from breathable and moisture transporting weaving fabric.

Indications:

- Shoulder dislocation.
- Post-surgical immobilization.
- Shoulder joint Pain.
- Forearm pain.
- Shoulder sprain.

Size	Inches	Cm
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-145



*Measure circumference around the chest below underarms



TD-303


Universal shoulder Immobilizer

Our Universal Shoulder Immobilizer comes in a flexible size to support shoulder immobilization.

Indications:

- Shoulder dislocation.
- Post-surgical immobilization.
- Shoulder joint Pain.
- Forearm pain.
- Shoulder sprain.

Size	Inches	Cm
CH	Below 26	Below 65
Universal	26-40	65-100
Special	40-52	100-130



*Measure circumference around the chest below underarms



TD-304


Clavicle Brace

Eight-shaped band design of Clavicle Brace is cushioned from soft polyurethane foam padding.

Indications:

- Clavicle fracture.
- Clavicle strain.
- Clavicle trauma.
- Poor posture problems.

Size	Inches	Cm
CH	22-28	56-71
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	132-147



*Measure circumference around the chest below underarms



TD-307


Clavicle Brace with Buckle

Eight-shaped Band Design of Clavicle Brace is Cushioned From Moisture Wicking Foam Fabric Padding.

Indications:

- Clavicle fracture.
- Clavicle strain.
- Clavicle trauma.
- Poor posture problems.

Size	Inches	Cm
CH	22-28	56-71
S	28-34	70-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	132-147



*Measure circumference around the chest below underarms



ARM & ELBOW SUPPORT



TD-301


Pouch Arm Sling (Tropical)

Our Pouch Arm Sling provides comfortable support to the wrist and arm with an adjustable strap with a buckle for proper immobilization.

Indications:

- Arm fracture.
- Wrist fracture.
- Arm & Wrist sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast.

Size	Inches	Cm
CH	22-28	56-71
S	28-34	71-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-147



*Measure circumference around the chest below underarms



TD-306


Pouch Arm Sling (Baggy)

Our baggy Pouch Arm Sling provides light-weight support to the wrist and arm through a thumb cradle, and an adjustable strap with a buckle for proper immobilization.

Indications:

- Arm Fracture.
- Wrist Fracture.
- Arm & Wrist Sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast rehabilitation.

Size	Inches	Cm
CH	22-28	56-71
S	28-34	71-85
M	34-40	85-100
L	40-46	100-115
XL	46-52	115-130
XXL	52-58	130-147



*Measure circumference around the chest below underarms



TD-504

Elbow Support

Our Elbow Support is made from a soft, elastic material. It provides compression and support for injured tendons.

Indications:

- Forearm/Elbow pain.
- Post Elbow
- Arthroscopy.
- Sports and workout stress.
- Post-cast rehabilitation.
- Postoperative weakness.
- Arthritis / Bursitis.

Size	Inches	Cm
S	7-10	18-25
M	10-12	25-30
L	12-14	30-35
XL	14-16	35-40



*Measure circumference around the elbow in extended position



TD-503

Tennis Elbow Support

Our Tennis Elbow support has adjustable straps and a reverse fastener for a strong grip and custom compression.

Indications:

- Tennis/Golf Elbow.
- Tendonitis.
- Swelling.

Size	Inches	Cm
S	7-10	18-25
M	10-12	25-30
L	12-14	30-35
XL	14-16	35-40



*Measure circumference around the elbow in extended position

ARM & ELBOW SUPPORT



TD-716


Arm Immobilizer **NEW/ARRIVAL**

An arm immobilizer is a useful medical device that can be used to support the healing of various arm injuries or surgeries. It can help to reduce pain and swelling. Prevent further injury, restrict the movement off the arm.


Indications:

- Fractures, Dislocations, and Soft Tissues injuries to Muscles, Ligaments, or Tendons
- Post-operative recovery of Arm
- Restrict the movement off Arm
- Immobilize the Elbow joint

Size	Inches	Cm
Universal	12.8-7.87	20.0-32.5



*Measure circumference around the elbow in extended position



TD-507

Forearm Splint

Our Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb.

Indications:

- Post-cast/Post-surgery protection and rehabilitation.
- Conservative management of Kienbock's disease.
- Rheumatic disease of the wrist.
- Arthritic hand positioning.

Size	Inches	Cm
Universal	7.5-10.50	19-26



*Measure distance from wrist to elbow



WRIST & FOREARM SUPPORT



TD-501

Wrist Brace

Designed with adjustable straps, moisture-wicking fabric and durable nylon reinforcement.

Indications:

- Sprains and strains.
- Mild Osteoarthritis.
- Swelling & inflammation.
- Prophylaxis for exercise and sports.

Size	Inches	Cm
S	4.9-6.1	12.5-15.5
M	6.1-7.3	15.5-18.5
L	7.3-8.5	18.5-21.5



*Measure circumference evenly around wrist joint



TD-502

Wrist Brace with Thumb Support

Designed with adjustable wrap-around straps, moisture-wicking fabric and durable nylon reinforcement.

Indications:

- Sprains and strains .
- Mild Osteoarthritis.
- Swelling & inflammation.
- Tendonitis.

Size	Inches	Cm
Universal	5.0-9.5	12.5-23.80



*Measure circumference evenly around wrist joint



TD-505


Wrist and Forearm Splint (Right/Left)

Our Wrist & Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb.

Indications:

- Post-cast management.
- Post-surgical protection and rehabilitation.
- Conservative management of Kienbock's disease.
- Rheumatic disease of the wrist.
- Arthritic hand positioning.

Size	Inches	Cm
S	7.2-8.5	18.3-21.6
M	8.5-9.8	21.6-24.9
L	9.8-11.1	24.9-28.2
XL	11.1-12.4	28.2-31.5



*Measure the circumference at a distance 6" from the wrist along the arm



TD-506

Elastic Wrist Splint (Right/Left)

Our athletic-inspired elastic splint is a lightweight wrist support with double straps and anatomical design for rigid immobilization, while maintaining the wrist in a functional position.

Indications:

- Carpal Tunnel Syndrome.
- Post-surgery.
- Post-cast management.
- Wrist sprains and strains.
- Tendonitis.

Size	Inches	Cm
S	4.9-6.1	12.5-15.5
M	6.1-7.3	15.5-18.5
L	7.3-8.5	18.5-21.5
XL	8.5-9.7	21.5-24.5



*Measure the circumference at a distance 6" from the wrist along the arm



FINGER SPLINTS



TD-601

Finger Cot

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort.

Indications:

- Finger injury.
- Post finger Fracture.
- Postoperative care.
- Sprains and bruises of fingers.

Size	Inches	Cm
S	< 2	< 5
M	2-3	5-7.5
L	3-4	7.5-10



*Length of fingers



TD-604

Finger Extension Splint

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort.

Indications:

- Metacarpal injury.
- Fracture.
- Post-surgery.
- Sprains and bruises of fingers.

Size	Inches	Cm
S	< 2	< 5
M	2-3	5-7.5
L	3-4	7.5-10



*Length of fingers



TD-602

Frog Splint

Made of light-weight but rigid aluminum frame to rehabilitate and align the post ligamentary damage of the Interphalangeal joints.


Indications:

- Interphalangeal (IP) joint Injury.
- Finger Injury.
- Finger sprain.
- Postoperative care.

Size	Inches	Cm
S	< 2	< 5
M	2-3	5-7.5
L	3-4	7.5-10



*Length of fingers



TD-605

Mallet Finger Splint

Ethafoam liner for cushioning and also suitable for post-operative care.

Indications:

- Bruising of the dorsum of the DIPJ.
- Mallet finger deformity.
- May help manage the extensor tendon injury.
- Minor fractures of the distal phalanges.
- Distal phalange postfracture care.
- Distal interphalangeal (DIP) joint conditions.
- Swollen fingertip.
- Post-surgical care.

Size	Inches	Cm
Universal	< 2-4	< 5-10



*Length of fingers

FINGER SPLINTS



TD-603

Thumb Spica Splint

A flexible and comfortable splint made of malleable material for easy application and removal.

Indications:

- Bursitis.
- Sprains of thumb.
- Avulsion fractures of the thumb.
- Soft tissue injuries.
- Thumb pain symptoms.
- CMC & MCP joint conditions.

Size	Inches	Cm
Universal	5-9.5	12.5-23.8



*Measure circumference evenly around wrist joint



TD-709

Garter Splint

NEW ARRIVAL

Garter splint is more adjustable than complete tights, allowing for a better fit and accommodating various fingers. It is more comfortable than traditional tights, especially for those with sensitive skin or those who experience discomfort from synthetic materials.

Indication

- Support and immobilization
- Control of pain
- Postoperative immobilization
- Prevention of further injury

Size	Inches	Cm
Universal	4.0-5.0	10-12.5

*Measure circumference evenly around wrist joint



TRACTION KITS



TD-801


Pelvic Binder

Our Pelvic Binder is made with ruggedness for extended and repetitive uses.

Indications:

- Lumbar spondylitis.
- Lower back pain.
- Postural fatigue.
- Postural deformities.
- Vertebral fractures.
- Mechanical stress due to poor posture.

Size	Inches	Cm
S	32-38	80-95
M	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3" below navel



TD-803


Pelvic Traction Kit (with Weight Bag)

A versatile traction kit that can be used in all kinds of tractions and on all kinds of beds.

Indications:

- Lumbar spondylitis.
- Lower back pain.
- Postural fatigue.
- Postural deformities.
- Vertebral fractures.

Size	Inches	Cm
S	32-38	80-90
M	38-44	95-110
L	44-50	110-125
XL	50-56	125-140



*Measure circumference 3" below navel



TD-803W

Traction Weight Bag

The Weight Bag applies uniform and controlled traction to the lumbar sacral, and thoracic areas of the vertebral column.

Indications:

- Convenient home use
- Controlled continuous traction
- Easy removal and application

Weight Capacity	
Liquid	10 Ltr.
Solid	12 Kg.



TD-804


Leg Traction Brace

Four sets of nylon based buckle closure with an additional locking mechanism to ensure a secure fit.

Indications:

- Hip/Sacral pain.
- Shaft fractures of the femur.
- Traction to relieve pain in the lower back
- When pelvic traction is not possible or not desired.

Size	Inches	Cm
S	8.5-11.0	21.50-27.90
M	11.0-13.50	27.90-34.20
L	13.5-16.0	34.20-40.60
XL	16.0-18.50	40.60-46.90



*Measure circumference around calf approx 6" below knee

TRACTION KITS



TD-802

Cervical Traction Kit + Weight Bag (Sleeping)

The Cervical Traction Kit (sleeping Type) is designed to provide traction to the cervical and upper dorsal vertebrae, relieving pressure by stretching the musculature and connective tissue components of the cervical spine. It expands the intervertebral joint space and can be used in a sleeping posture, providing the convenience of traction therapy at home.

Indications:

- Cervical injuries
- Cervical disc herniation
- Cervical radiculopathy
- Cervical muscle strain

Size	Universal
------	-----------



TD-719

Cervical Traction Kit + Weight Bag (Sitting)

The Cervical Traction Kit (Sitting Type) is used to relieve pressure on muscle, nerve tissue, and minor fractures of the cervical spine. It works by stretching the cervical spine's musculature and connective tissue components, expanding the intervertebral joint space, and reducing pressure on the cervical spine.

Indications:

- Cervical Disc Herniation
- Cervical Radiculopathy
- Cervical Muscle Strain
- Neck Pain And Stiffness

Size	Universal
------	-----------



TD-805

Skin Traction Set (Puf Liner)

Our Skin Traction Kit is Made With Foam Liner Designed To Ensure There is No Friction Damage To The Skin And No Slippage During Traction.

Indications:

- Need For Indirect Traction Tto Extremities
- Pain & Spasms

Size	Inches	Cm
Child	< 12	< 30
Adult	> 12	> 30



THIGH, KNEE, ANKLE & CALF SUPPORT



TD-416

Thigh Support

The weaving method for our Thigh Support provides seamless stretching stockings with a tapered top to bottom design to hold in a natural shape.

Indications:

- Hamstring strain.
- Muscle spasms or cramping.
- Arthritis.
- Stiffness or trauma.
- Pain & inflammation.
- Sports Injury.

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6-22	49-55
XL	22-24.4	55-61



*Measure circumference around mid thigh approx 6" above knee



TD-403


Knee Cap

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6-22.0	49-55
XL	22.0-24.4	55-61
XXL	24.4-26.8	61-67
XXXL	26.8-29.2	67-73



*Measure circumference around mid thigh approx 6" above knee



TD-404


Knee Cap Open Patella

Spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee



TD-405


Knee Cap Rigid Hinge

An anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee Sprain.
- Weak knees.

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee

THIGH, KNEE, ANKLE & CALF SUPPORT



Elastic Knee Support


Made from elastic, comfortable and durable material with adjustable strap for proper compression and customized fit.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee Sprain.
- Weak knees.

TD-409

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee



Patellar Support

Provides enhanced kinesthesia for promoting blood circulation and healing with customized compression and padding.

Indications:

- Patellofemoral pain syndrome.
- Patellar tendonitis and arthritis.
- Osgood-Schlatter disease.
- Patellar dislocation and subluxation.
- Mild instabilities of the knee joint.

TD-414

Size	Inches	Cm
Universal	12-20	30-50



*Measure circumference around the knee



Knee Immobilizer (14")


Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

Indications:

- Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- Acute sprains/strains of the knee.
- Post-cast care.
- Post/Preoperative care.

TD-406

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee



Knee Immobilizer (19")


Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

Indications:

- Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- Acute sprains/strains of the knee.
- Post-cast care.
- Post/Preoperative care.

TD-407

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee

THIGH, KNEE, ANKLE & CALF SUPPORT



TD-408


Knee Immobilizer (22")

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

Indications:

- Ligament injuries.
- Tibial plateau fractures.
- Osteochondral injury.
- Patellar tendon injury.
- Acute sprains/strains of the knee.
- Post-cast care.
- Post/Preoperative care.

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee



TD-411

R.O.M Knee Brace (22")

Our R.O.M knee brace provides locking at any angle or motion between any two angles with intervals of 15 degrees.

Indications:

- Arthritis/Osteoarthritis.
- Hamstring Injury.
- Knee ligament instability.
- Knee sprain.
- Lateral Collateral Swelling & inflammation.
- Total knee replacement.

Size	Inches	Cm
Universal	14.8-22	37-55



*Measure circumference around mid thigh approx 6" above knee



TD-417


Calf Support

Hypoallergenic elastic fabric from unique fibres blended with cotton prevents skin irritations and perfect compression.

Indications:

- Weak or injured calves, fatigued legs, and stressed muscles.
- Aching calf muscles.
- Cramps and painful calf conditions.
- Muscle strain.

Size	Inches	Cm
S	9.8-11.3	25-29
M	11.3-12.8	29-33
L	12.8-14.3	33-36
XL	14.3-15.8	36-40



*Measure circumference approx 2" below knee



TD-413


Foot Drop Splint (Right/Left)

Useful for patients under critical care of foot injury. Lightweight, robust and flexible polypropylene shell in one piece.

Indications:

- Muscle weakness due to stroke, spinal cord injury, muscular dystrophy, cerebral paralysis, poliomyelitis, multiple sclerosis and peripheral neuropathy.
- Recovery from surgical or neurological treatments for foot drop.
- Muscle tightening and runner's cramp.
- Plantar fasciitis.
- Achilles tendonitis.

Size	Europe	USA	MM
CH	25-30	8-12	157-190
S	34-36	2.75-4.35	217-230
M	37-39	5.25-6.75	237-250
L	40-42	7.50-9.00	257-270
XL	43-45	9.75-11.5	277-291



*Measure Shoe Size

THIGH, KNEE, ANKLE & CALF SUPPORT



Ankle Br ace

Easy to wear and remove brace made of rigid polypropylene material for strong immobilization.

Indications:

- Mild ankle sprains
- Recurring ankle injury.
- Tarsal tunnel syndrome.
- Swollen ankle.
(Not indicated for severe ankle instabilities)

TD-412

Size	Inches	Cm
CH	Below 7.2	Below 18
S	7.2-8.4	18-21
M	8.4-9.6	21-24
L	9.6-10.8	24-27
XL	10.8-12	27-30



*Measure circumference approx 2" above the ankle joint



Ankle Binder

The figure of eight design makes it comfortably wearable, even with shoes.

Indications:

- Supination injuries.
- Post-supination injuries.
- Syndesmotoc ligament injury.
- Chronic ligament instability in the upper & lower ankle.

TD-401

Size	Inches	Cm
S	7.2-8.4	18-21
M	8.4-9.6	21-24
L	9.6-10.8	24-27
XL	10.8-12	27-30



*Measure circumference approx 2" above the ankle joint



Anklet

Provides compression for the foot and ankle joint.

Indications:

- Lymphoedema.
- Lipedema.
- Chronic venous insufficiency.
- Varicose veins.
- Post Thrombotic Syndrome.
- Venous Stasis Ulcer.
- Post Sclerotherapy.
- Lipodermatosclerosis.

TD-402

Size	Inches	Cm
S	7.2-8.4	18-21
M	8.4-9.6	21-24
L	9.6-10.8	24-27
XL	10.8-12	27-30



*Measure circumference approx 2" above the ankle joint



COMPRESSION GARMENTS



Compression Stocking Mid Thigh


Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

TD-703

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference a mid thigh approx 6" above knee



Compression Stocking Below Knee


Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort.

Indications:

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

TD-704

Size	Inches	Cm
S	8.5-11.0	21.5-27.90
M	11.0-13.50	29.90-34.20
L	13.50-16.0	34.20-40.60
XL	16.0-18.50	40.60-46.90
XXL	18.50-21.0	46.90-52.50



*Measure circumference around calf approx 6" below knee



Medical Compression Stocking Mid Thigh-Classical

Medical compression stockings are specially designed to apply pressure to your legs, helping to maintain blood flow and reduce discomfort and swelling.

Indications:

- Maintain blood flow
- Reduce discomfort & swelling
- Improving lymphatic drainage
- Managing orthostatic hypotension

TD-710

Perimeter/Model	S	M	L	XL
cB	20-22	23-25	26-29	30-32
cC	30-39	33-41	36-45	39-48
cD	29-38	32-40	35-44	38-47
cE	49-63	54-66	58-71	60-75



*Measure circumference a mid thigh approx 6" above knee



Medical Compression Stocking Below Knee-Classical


Medical compression stockings are specially designed to apply pressure to your legs, helping to maintain blood flow and reduce discomfort and swelling.

Indications:

- Maintain blood flow
- Reduce discomfort & swelling
- Improving lymphatic drainage
- Managing orthostatic hypotension

TD-711

Perimeter/Model	S	M	L	XL
cB	20-22	23-25	26-29	30-32
cC	30-39	33-41	36-45	39-48
cD	29-38	32-40	35-44	38-47
cE	49-63	54-66	58-71	60-75



*Measure circumference around calf approx 6" below knee

COMPRESSION GARMENTS



TD-714

Prosthetic Compression Stocking - Mid Thigh

Prosthetic compression stockings prevent edema, leg discomfort, tired aching, deep vein thrombosis, and mild swelling. They work by applying pressure to the legs, which helps to improve blood flow and prevent blood from pooling and clotting.

Indications:

- Prevention of deep vein thrombosis
- Venous disorders such as edema, phlebitis, and thrombosis
- Reduction of pain and swelling in the legs
- Improvement of residual limb blood circulation




TD-715

Prosthetic Compression Stocking - Below Knee

Prosthetic compression stockings prevent edema, leg discomfort, tired aching, deep vein thrombosis, and mild swelling. They work by applying pressure to the legs, which helps to improve blood flow and prevent blood from pooling and clotting.

Indications:

- Prevention of deep vein thrombosis
- Venous disorders such as edema, phlebitis, and thrombosis
- Reduction of pain and swelling in the legs
- Improvement of residual limb blood circulation



Above knee lengths:


- Extra Long 14" (35 cm)
- Long 12" (30 cm)
- Medium 10" (25 cm)
- Short 8" (20 cm)

Size (cm, Inches)

	XS	S	M	L	XL	2XL	3XL	4XL	5XL	6XL	7XL
cB	19-21 cm 7 1/2 - 8 1/4	21-23 cm 8 1/4 - 9	23-25 cm 9 - 9 3/4	25-27 cm 9 3/4 - 10 3/4	27-29 cm 10 3/4 - 11 1/2	28-30 cm 11 - 11 1/2	30-32 cm 11 1/2 - 12 1/4	32-35 cm 12 1/4 - 13 3/4	35-38 cm 13 3/4 - 15	37-40 cm 14 1/4 - 16	39-42 cm 15 - 16 1/2
cC	29-32 cm 11 1/2 - 12 1/4	32-35 cm 12 1/4 - 13 3/4	35-38 cm 13 3/4 - 15	39-42 cm 15 1/4 - 16 1/4	42-45 cm 16 1/4 - 17 1/4	44-47 cm 17 1/4 - 18 1/4	46-49 cm 18 - 19 1/4	49-53 cm 19 1/4 - 20 3/4	52-56 cm 20 3/4 - 22	56-62 cm 22 - 22 3/4	60-66 cm 23 3/4 - 26
cD	28-31 cm 11 - 12 1/4	31-34 cm 12 1/4 - 13 1/4	34-37 cm 13 1/4 - 14 1/4	38-41 cm 15 - 16 1/4	41-44 cm 16 1/4 - 17 1/4	43-46 cm 17 - 18	45-48 cm 17 1/4 - 18 1/4	48-52 cm 18 1/4 - 20 3/4	51-55 cm 20 - 21 1/4	55-61 cm 21 1/4 - 24	59-65 cm 23 1/4 - 25 3/4
cE	41-44 cm 16 1/4 - 17 1/4	43-47 cm 17 - 18 1/4	46-50 cm 18 - 19 1/4	49-53 cm 19 1/4 - 20 3/4	52-56 cm 20 3/4 - 22	55-59 cm 21 1/4 - 23 1/4	58-62 cm 22 1/4 - 24 1/4	62-67 cm 24 1/4 - 26 1/4	65-71 cm 25 3/4 - 28	70-77 cm 27 1/4 - 30 1/4	76-85 cm 30 - 33 3/4
cG	42-48 cm 16 1/4 - 19 1/4	46-52 cm 18 - 20 3/4	51-57 cm 20 - 22 3/4	56-62 cm 22 - 24 3/4	60-66 cm 23 3/4 - 26	64-70 cm 25 1/4 - 27 1/4	68-74 cm 26 1/4 - 29 1/4	73-80 cm 28 1/4 - 31 1/4	77-84 cm 30 1/4 - 33	82-89 cm 32 1/4 - 35	87-95 cm 34 1/4 - 37 1/4
A-ID	34-36 cm 13 1/4 - 14 1/4	36-38 cm 14 1/4 - 15	37-39 cm 14 3/4 - 15 1/4	39-41 cm 15 1/4 - 16 1/4	40-42 cm 15 3/4 - 16 3/4	42-44 cm 16 3/4 - 17 1/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4
A-IG	63-67 cm 24 3/4 - 26 1/4	66-70 cm 26 - 27 1/4	69-73 cm 27 1/4 - 28 3/4	72-76 cm 28 1/4 - 29	75-79 cm 29 1/4 - 31 1/4	78-82 cm 30 3/4 - 32 1/4	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33

EU Standard

■ Class1 (18-21mmHg) ■ Class2 (23-32mmHg) ■ Class3 (34-46mmHg)



Below knee lengths:

- Extra Long 22" (55 cm)
- Long 20" (50 cm)
- Medium 18" (45 cm)
- Short 16" (40 cm)

Size (cm, Inches)

	XS	S	M	L	XL	2XL	3XL	4XL	5XL	6XL	7XL
cB	19-21 cm 7 1/2 - 8 1/4	21-23 cm 8 1/4 - 9	23-25 cm 9 - 9 3/4	25-27 cm 9 3/4 - 10 3/4	27-29 cm 10 3/4 - 11 1/2	28-30 cm 11 - 11 1/2	30-32 cm 11 1/2 - 12 1/4	32-35 cm 12 1/4 - 13 3/4	35-38 cm 13 3/4 - 15	37-40 cm 14 1/4 - 16	39-42 cm 15 - 16 1/2
cC	29-32 cm 11 1/2 - 12 1/4	32-35 cm 12 1/4 - 13 3/4	35-38 cm 13 3/4 - 15	39-42 cm 15 1/4 - 16 1/4	42-45 cm 16 1/4 - 17 1/4	44-47 cm 17 1/4 - 18 1/4	46-49 cm 18 - 19 1/4	49-53 cm 19 1/4 - 20 3/4	52-56 cm 20 3/4 - 22	56-62 cm 22 - 22 3/4	60-66 cm 23 3/4 - 26
cD	28-31 cm 11 - 12 1/4	31-34 cm 12 1/4 - 13 1/4	34-37 cm 13 1/4 - 14 1/4	38-41 cm 15 - 16 1/4	41-44 cm 16 1/4 - 17 1/4	43-46 cm 17 - 18	45-48 cm 17 1/4 - 18 1/4	48-52 cm 18 1/4 - 20 3/4	51-55 cm 20 - 21 1/4	55-61 cm 21 1/4 - 24	59-65 cm 23 1/4 - 25 3/4
cE	41-44 cm 16 1/4 - 17 1/4	43-47 cm 17 - 18 1/4	46-50 cm 18 - 19 1/4	49-53 cm 19 1/4 - 20 3/4	52-56 cm 20 3/4 - 22	55-59 cm 21 1/4 - 23 1/4	58-62 cm 22 1/4 - 24 1/4	62-67 cm 24 1/4 - 26 1/4	65-71 cm 25 3/4 - 28	70-77 cm 27 1/4 - 30 1/4	76-85 cm 30 - 33 3/4
cG	42-48 cm 16 1/4 - 19 1/4	46-52 cm 18 - 20 3/4	51-57 cm 20 - 22 3/4	56-62 cm 22 - 24 3/4	60-66 cm 23 3/4 - 26	64-70 cm 25 1/4 - 27 1/4	68-74 cm 26 1/4 - 29 1/4	73-80 cm 28 1/4 - 31 1/4	77-84 cm 30 1/4 - 33	82-89 cm 32 1/4 - 35	87-95 cm 34 1/4 - 37 1/4
A-ID	34-36 cm 13 1/4 - 14 1/4	36-38 cm 14 1/4 - 15	37-39 cm 14 3/4 - 15 1/4	39-41 cm 15 1/4 - 16 1/4	40-42 cm 15 3/4 - 16 3/4	42-44 cm 16 3/4 - 17 1/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4	43-45 cm 17 - 17 3/4
A-IG	63-67 cm 24 3/4 - 26 1/4	66-70 cm 26 - 27 1/4	69-73 cm 27 1/4 - 28 3/4	72-76 cm 28 1/4 - 29	75-79 cm 29 1/4 - 31 1/4	78-82 cm 30 3/4 - 32 1/4	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33	80-84 cm 31 1/4 - 33

EU Standard

■ Class1 (18-21mmHg) ■ Class2 (23-32mmHg) ■ Class3 (34-46mmHg)



FOOT, HEEL & TOES SUPPORT



TD-706


Arch Support

Silicone based material absorbs shock, reduces friction and improves general foot comfort. Controls foot position and prevents arch pain.


Indications:

- Flat foot and associated problems.
- Walking discomfort.
- Plantar fasciitis.

Size	European	American
S	29-34	CH 11-AD 3
M	34-39	3.0-7.0
L	39-43	7.0-10.50



*Measure Shoe Size



TD-705

Toe Separator

Odour-free, silicone based toe separator helps release pressure and realign the toes.

Indications:

- Bunion, corns and calluses pain.
- Hallux Valgus.
- Toe misalignment.
- Plantar Fasciitis.
- Crooked toes, claw toes & hammertoes.

Size	European	American
S	34-36	2.8-4.4
M	37-39	5.3-6.8
L	40-42	7.5-9.0



*Measure Shoe Size



TD-702

Insole

Made of medical grade silicone for durability and comfort.


Indications:

- Achilles Tendonitis.
- Heel pain.
- Heel spurs.
- Pain due to walking.

Size	European	American
S	34-36	2.8-4.4
M	37-39	5.3-6.8
L	40-42	7.5-9.0
XL	43-45	9.8-11.5



*Measure Shoe Size



TD-701

Heel Cushion

Made of medical grade silicone for durability and comfort.

Indications:

- Achilles Tendonitis.
- Heel pain.
- Heel spurs.

Size	European	American
S	34-36	2.8-4.4
M	37-39	5.3-6.8
L	40-42	7.5-9.0
XL	43-45	9.8-11.5



*Measure Shoe Size

FOOT, HEEL & TOES SUPPORT



TD-707

Metatarsal Pad

Decreasing the pressure on the ball of your foot by supporting the metatarsal bone.

Indications:

- Plantar fasciitis
- Plantar fat pad atrophy
- Metatarsalgia (stone bruise)

Size	European	American
Universal	34-45	2.8-11.75

*Measure Shoe Size



TD-418


Cast Shoe

A light weight cast shoe which improves balance and helps early ambulation.

Indications:

- Neuropathic diabetic foot ulcers
- Forefoot and mid-foot injury
- Post ligament, soft tissues and tendon surgery

Size	European	American
S	34-36	2.8-4.4
M	37-39	5.3-6.8
L	40-42	7.5-9.0
XL	43-45	9.8-11.5



*Measure Shoe Size



TD-708

Bunion Splint

NEW ARRIVAL

The anatomically moulded splint is cast from high-quality polymer for rigidity and durability. Fastened metal hinges for improved durability. Polyethylene foam provides good cushioning and added comfort washable, ventilated and adapts to nearly any size to easily wear with shoes.

Indication

- Foam padding offers enhance comfort
- Adjustable shape
- Easy to clean
- Can be worn day or night

Size	European	American	MM
Universal	34-35	2.8-11.75	217-291

NEOPRENE PRODUCTS



TD-305


Shoulder Support (Neoprene)

Made from Neoprene bonded nylon fabric which is four-way stretchable with an anatomical shape to ensure perfect fitting and optimal immobilization.

Indications:

- Shoulder Injury.
- Rotator cuff injury.
- Shoulder dislocation
- Shoulder inflammation.
- Arthritic pain of the shoulder joint.
- Bursitis and Tendonitis.

Size	Inches	Cm
Universal	26-40	65-100
Special	40-52	100-130



*Measure circumference around the chest below underarms



TD-112

Waist Belt (Neoprene)

Our wrap-around design waist belt is made with neoprene bonded fabric for four-way stretchability and maximum support.

Indications:

- Weak abdominal muscles.
- Swelling & stiffness in abdominal muscles.
- Postoperative, post-delivery support.

Size	Inches	Cm
Universal	28-50	70-125



*Measure circumference around the navel



TD-508

Wrist Brace (Neoprene)

Created out of triple-layered neoprene. The top layer features continuous loops with rugged nylon fabric in the inner layer.

Indications:

- Tendonitis.
- Tenosynovitis.
- Arthritic care of wrist.
- Sprains and strains.
- Dislocation of the wrist joint.
- Postoperative and postcast care.
- Prophylaxis for exercise and sports.

Size	Inches	Cm
Universal	5.0-9.5	12.5-23.80



*Measure circumference evenly around wrist joint



TD-509


Wrist Brace With Thumb Support (Neoprene)

Neoprene based design with adjustable wraparound straps, moisture-wicking fabric and durable nylon reinforcement.

Indications:

- Ligamentous laxity of the thumb.
- Median nerve injury.
- Skier's thumb.
- Carpometacarpal joint injury.
- Post operative and post-cast care.
- Prophylaxis for exercise and sports.

Size	Inches	Cm
Universal	5.0-9.5	12.5-23.80



*Measure circumference evenly around wrist joint

NEOPRENE PRODUCTS



Knee Wrap (Neoprene)

Made with triple-layer neoprene fabric. The triple hook loop closing system ensures easy wearing, removal and proper fit.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees.

TD-415



Knee Wrap Hinged (Neoprene)


Simple and quick put on with an easy wrap-around application made from tapered, four-way stretchable neoprene material.

Indications:

- Arthritic pain.
- Osteoarthritis of the knee.
- Ligaments or tendon injury.
- Quad cramps.
- Knee sprain.
- Weak knees.


TD-410

Size	Inches	Cm
Universal	14.8-22	37-55
Special	22-29.2	55-73



*Measure circumference around mid thigh approx 6" above knee

Size	Inches	Cm
S	14.8-17.2	37-43
M	17.2-19.6	43-49
L	19.6 - 22	49-55
XL	22-24.4	55-61
XXL	24.4-26.8	61-67



*Measure circumference around mid thigh approx 6" above knee

PHYSIO AIDS



Exercise Band (Actiband)

It can be used for a variety of applications. Used particularly lower body exercise to increase strength, body Balance and Strengthen muscles.

Indications:

- Improve strength & flexibility in glutes, Knees and hips.
- Convenient and Portable.
- Powder-free surface for better gripping.

Resistance Level	Colour Code	Resistance in Lbs*	Resistance in Kgs*
Light	Yellow	09.98	4.5
Medium	Red	11.86	5.4
Strong	Green	13.73	6.2
Extra Strong	Blue	15.60	7.0
Super Strong	Black	16.85	7.7

TD-1301






Exercise Gel Ball

Anti-stress exercise gel ball in three colours depending on hardness. Good for building strength of fingers, palm and wrist.

Indications:

- Post fracture.
- Stress.
- Weak forearm and grip strength.

Soft	Medium	Hard
		

TD-902



Weight Cuff

Our weight cuff is fit for graded exercise therapy (GET) and also cosmetic arm muscle shaping.

Indications:

- Graded exercise therapy (GET).
- Poor muscle tone.
- Post-cast care.

Size	½ Kg	1 Kg	2 Kg
Variations of Weight Cuff fit all wrists and ankles.			

TD-901



WELLNESS PRODUCTS



Back Rest Back Rest

Promotes Good Posture And Correct Spinal Alignment By Allowing The Spine To Follow Natural Curve.

Indications:

- Postoperative discomfort.
- Post-injury pain.
- Poor posture.
- Lower back pain.
- Pregnancy discomfort.

TD-1101

Size	Cm
Universal	152-192



Back Rest -Lumbo **NEW ARRIVAL**

The purpose of the Tesca backrest is to provide support and comfort in various settings, such as beds, chairs, and office furniture. This backrest contributes to good posture, supporting the natural curve of the spine and reducing the risk of back pain associated with prolonged sitting. It is a well-designed backrest that promotes ergonomic sitting positions, enhances comfort, ensures proper alignment and reduces strain on the back, neck, and shoulders.

Indication

- Support and Stability
- Posture Improvement
- Prevention of Musculoskeletal Problems
- Enhanced Comfort

TD-729

Size	Inches (LXBXH)	Cm (LXBXH)
Universal	16.4X16.5X2.7	41.7X42X7



Coccyx Cushion Seat

The high-density foam of our coccyx cushion features a U-shaped cut that gently suspends your coccyx to relieve pressure on your spinal cord.

Indications:

- Pudendal nerve entrapment pain.
- Chronic pelvic pain syndrome.
- Rectal surgery.
- Prostatitis.
- Anal surgery.
- Coccyx fracture.
- Hemorrhoids.
- Interstitial Cystitis.
- Pelvic floor pain.
- Pressure sores.

TD-1102

Size	Inches (LXBXH)	Cm (LXBXH)
Universal	16.4X16.5X2.7	41.7X42X7



Orthopaedic Donut Pillow **NEW ARRIVAL**

Orthopedic donut pillow can alleviate direct pressure and provide high-quality support. They can adapt to your body shape to work for you and provide effective pain relief. It is not just for those suffering from problems associated with sitting. They are the perfect solution for those who have injuries or discomfort in the lower region, for those who suffer from hemorrhoids or even for women after giving birth.

Indication

- Relief from Pain and Discomfort
- Improve Posture
- Pressure Relief
- Support For Sitting
- Promotion of Spine Alignment

TD-717

Size
Universal

WELLNESS PRODUCTS



Travel Neck Pillow **NEW ARRIVAL**

Travel neck pillow provide support and maintain the correct posture, reducing neck strain and preventing fatigue. It align your body to the seat, promoting better posture and reducing the risk of muscle cramps. It also help align your head and neck, allowing to open up and reducing snoring.

size	Neck Height Inches	Neck Height Cm
Uni.	2.5 - 5.0	6.25 - 12.5

Indication

- Support and Comfort
- Posture Alignment
- Prevention of Tiredness
- Prevention of Neck Pain

TD-718



Tourniquet

Our tourniquet is padded with soft foam on The inner surface and constructed with A lasting hook & loop strap with a tough Polypropylene fastener.

Size	Inches	Cm
Universal	9.8-11.10	24.90-28.20




*Measure circumference around the elbow in extended position

Indications:

- Stop constant bleeding.
- Restrain the circulation of blood to an extremity for a period of time.

TD-1201



Orthopaedic Heating Pad **NEW ARRIVAL**

An orthopaedic heating pad is designed to provide relief from musculoskeletal pain and discomfort. It is often used for conditions such as back pain, arthritis, and muscle aches. These heating pads typically feature au to heat control technology, fast heating, and ultra-soft materials for comfort.

Size	Cms
Regular	29.5X21.5
Extra-Large	34X25

Indication

- 3 Comfortable Heat Levels
- Made Of Soft Materials
- Low Power Consumption
- Temperature Control

TD-720



HOT & COLD THERAPY AIDS



TD-1001

Hot & Cold Pack

Pack Made Of Heating Refrigerant Gel Which Instantly Relieves Muscles Pain And Other Aches Due To Wound, Sprain, Etc.

Indications:

- Lumbago (Low Back Pain).
- Joint pain & swelling.
- Muscle spasms.

Size

Universal



OTHER PRODUCTS



TD-1402/03/04/05

ACTICREPE (Cotton Crepe Bandage B.P.)

Tesca Crepebandages are to provide support, padding, and protection in various situations, such as sprains, strains, muscle injuries, and wounds. It is designed with high elasticity, allowing it to conform well to the body and provide light support while allowing movement. It helps to treat varicose sprains, muscle injuries, and fractures, offering compression and stability to the affected area.

Indication:

- Comfort and compression
- Injury support
- Emergency first aid dressing
- Reducing swelling effectively

Product Code	Product Size
TD-1402	6 Cm x 4 Mtr
TD-1403	8 Cm x 4 Mtr
TD-1404	10 Cm x 4 Mtr
TD-1405	15 Cm x 4 Mtr



TD-3209

Disposable 3 Ply Face Mask

3-ply mask is a disposable mask that has 3 layers of non-woven fabric to help you breathe clean air.

Indications:

- Adjustable nose bridge
- Seamless crimping & breathability

Size

Universal



TD-1409/10

ACTIPOP (Plaster of Paris Bandage)

“ACTIPOP” Plaster of Paris is Made From Specially Woven Leno Gauze Evenly Impregnated With Finest Quality Plaster of Paris Powder That Product An Even Cast And Excellent Finish. It is Used By Orthopedic Surgeons/Technicians for Immobilization of Simple & Compound Fractures.

Product Features:

- Fast Setting Time
- High Load Bearing Cast
- Green Setting Time
- No Dry Patches
- Minimum Powder Loss
- Fast Saturation Time
- Water Absorption Time Is 2 To 4 Sec
- Fully Creamy Bandage
- Made By Leno Gauze Net Cotton Yarn
- High Quality Gymsum

Product Code	Product Size
TD-1409	10cm X 2.7 mtr
TD-1410	15cm X 2.7 mtr




NOTES

A series of horizontal dashed lines for writing notes.



NOTES

A series of horizontal dashed lines for writing notes.



CONTACT SALES



Mob./WhatsApp: +91-9829132777
Email: www.tescaglobal.com
Website: www.tescaglobal.com

Address:
Tesca Technologies Pvt. Ltd.
IT-2013, Ramchandrapura Industrial Area,
Sitapura Extension, Jaipur-302022, India
Mob: +91-9829132777, 9413330765
Email: info@tesca.in



SCAN ME